

Anti-PGP9.5 (7B8) Mouse antibody



Catalog Number: 168282

Product name

Anti-PGP9.5 (7B8) Mouse antibody

Specificity

Human, Monkey

Antibody description

PGP9.5 (7B8) Mouse monoclonal antibody

Preparation

Antigen: Purified recombinant human UCHL1 / PGP9.5 protein fragments expressed in E.coli.

Formulation

Purified mouse monoclonal in buffer containing 0.1M Tris-Glycine (pH 7.4 150 mM NaCl) with 0.02% sodium azide 50% glycerol

Storage

Store at 4°C short term. Store at -20°C long term. Avoid freeze / thaw cycle.

Clonality

Monoclonal

Ig Type

Mouse IgG2b

Applications


WB, ICC/IF

Dilutions


WB: 1/1000

;ICC: 1/300

Validations



Immunocytochemistry stain of COS7 using UCHL1 / PGP9.5 mouse mAb (1:300).



Western blot detection of UCHL1 / PGP9.5 in U251,IMR32,Raji and COS7 cell lysates and using UCHL1 / PGP9.5 mouse mAb (1:1000 diluted).Predicted band size: 25KDa.Observed band size: 28KDa.