

# MLL5 antibody



Catalog Number: 112682

## Product name

MLL5 antibody

## Immunogen

[Human MLL5 Recombinant protein \(GST tag\)](#)

## Specificity

human, mouse; other species not tested.

## Antibody description

MLL5 Rabbit Polyclonal antibody. Positive WB detected in HL-60 cells, HeLa cells, Jurkat cells, mouse brain tissue. Observed molecular weight by Western-blot: 100 kDa, 70 kDa

## Preparation

This antibody was obtained by immunization of MLL5 recombinant protein (Accession Number: BC062583). Purification method: Antigen affinity purified.

## Formulation

PBS with 0.1% sodium azide and 50% glycerol pH 7.3.

## Storage

Store at -20°C. DO NOT ALIQUOT

## Clonality

Polyclonal

## Ig Type

Rabbit IgG


## Applications

ELISA, WB

## Dilutions

Recommended Dilution:  
WB: 1:200-1:2000

## Validations



HL-60 cells were subjected to SDS PAGE followed by western blot with Catalog No:112682(MLL5 antibody) at dilution of 1:800