

Human TAC1/NKA/NKNA/TAC2 (His tag) recombinant protein



Catalog Number: 519682

General Information

Gene Name Synonym

Protachykinin-1, PPT, TAC1

Protein Construction

A DNA sequence encoding the Human TAC1/NKA/NKNA/TAC2 Glu20-Arg129 is expressed with a 6His tag at the N-terminus.

Organism

Human

Expression Host

E. coli

QC Testing

Purity

>95% as determined by SDS-PAGE.

Endotoxin

<1.0 EU per μ g protein as determined by the LAL method.

Stability

Store for up to 6 months at -20°C

Formulation


Supplied as a 0.2 μ m filtered solution of 20mM Tris, 150mM NaCl, pH 8.0 (10% glycerol 2M Urea).

Usage Guide

Storage

Store at <-20°C, minimize freeze-thaw cycles.

SDS-PAGE



Human TAC1/NKA/NKNA/TAC2 (His tag)
recombinant protein