

# HDAC7 antibody

Catalog Number: 111382

## Product name

HDAC7 antibody

## Immunogen

[Human HDAC7 Recombinant protein \(His tag\)](#)

## Specificity

Human; other species not tested.

## Antibody description

HDAC7 Rabbit Polyclonal antibody. Positive IP detected in K-562 cells. Observed molecular weight by Western-blot: 102 kDa

## Preparation

This antibody was obtained by immunization of HDAC7 recombinant protein (Accession Number: BC006453). Purification method: Antigen affinity purified.

## Formulation

PBS with 0.1% sodium azide and 50% glycerol pH 7.3.

## Storage

Store at -20°C. DO NOT ALIQUOT

## Clonality

Polyclonal

## Ig Type

Rabbit IgG


## Applications

ELISA, IP

## Dilutions

Recommended Dilution:  
IP: 1:200-1:1000

## Validations



IP Result of anti-HDAC7 (IP:Catalog No:111382, 4ug; Detection:Catalog No:111382 1:300) with K-562 cells lysate 4000ug.