

# NUP62CL antibody



Catalog Number: 113382

## Product name

NUP62CL antibody

## Specificity

Human; other species not tested.

## Antibody description

NUP62CL Rabbit Polyclonal antibody. Positive WB detected in HeLa cells, Jurkat cells, SH-SY5Y cells. Observed molecular weight by Western-blot: 30kDa

## Preparation

This antibody was obtained by immunization of NUP62CL recombinant protein (Accession Number: BC017799). Purification method: Antigen affinity purified.

## Formulation

PBS with 0.1% sodium azide and 50% glycerol pH 7.3.

## Storage

Store at -20°C. DO NOT ALIQUOT

## Clonality

Polyclonal

## Ig Type

Rabbit IgG

## Applications


ELISA, WB

## Dilutions

Recommended Dilution:

WB: 1:200-1:2000

## Validations



HeLa cells were subjected to SDS PAGE followed by western blot with Catalog No:113382(NUP62CL Antibody) at dilution of 1:600