

# FSD2 antibody



Catalog Number: 110782

---

## Product name

FSD2 antibody

## Specificity

Human, Mouse, Rat; other species not tested.

## Antibody description

FSD2 Rabbit Polyclonal antibody. Positive WB detected in mouse heart tissue, rat heart tissue. Observed molecular weight by Western-blot: 95-105 kDa

## Preparation

This antibody was obtained by immunization of FSD2 recombinant protein (Accession Number: XM\_024449845). Purification method: Antigen affinity purified.

## Formulation

PBS with 0.1% sodium azide and 50% glycerol pH 7.3.

## Storage

Store at -20°C. DO NOT ALIQUOT

## Clonality

Polyclonal

## Ig Type

Rabbit IgG

## Applications


ELISA, WB

## Dilutions

Recommended Dilution:

WB: 1:200-1:2000

## Validations



mouse heart tissue were subjected to SDS PAGE followed by western blot with Catalog No:110782(FSD2 Antibody) at dilution of 1:600