

Anti-VIII antibody



Catalog Number: 175092

Product name

Mouse IgG1

Anti-VIII antibody

Applications

Specificity

ELISA

M13 phage

Dilutions

Antibody description

ELISA: 1:5,000-1:20,000

Mouse monoclonal antibody to VIII

Validations

Preparation

This antigen of this antibody was synthetic peptide.

Formulation

Liquid, 1*PBS (pH7.4).

Storage

Store at +4°C after thawing. Aliquot store at -20°C. Avoid repeated freeze / thaw cycles.

Clonality

Monoclonal

Ig Type

Elisa Test of Anti-M13-HRP antibody Binding to M13 Phage




Fig1.; Elisa Test of M13 antibody (HRP) Binding to M13 Phage.; Samples were blocked with 1% BSA for 1 hour at 37°C and incubated with the primary antibody (1/10000) for 1 hour at 25°C. The M13 antibody was repeatedly frozen and thawed for 4 times at -80°C, 4°C, 37°C.