

Anti-MRPL18 antibody



Catalog Number: 176092

Product name

Anti-MRPL18 antibody

Specificity

Human

Antibody description

Rabbit polyclonal antibody to MRPL18

Preparation

This antigen of this antibody was the antiserum was produced against synthesized peptide derived from human mrpl18. aa range:1-50

Formulation

Liquid, Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.

Storage

Store at +4°C after thawing. Aliquot store at -20°C or -80°C. Avoid repeated freeze / thaw cycles.

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications

WB, IHC-P

Dilutions

WB: 1:500-1:2000

IHC-P: 1:100-1:300

Validations



Fig1:; Western Blot analysis of various cells using MRP-L18 Polyclonal Antibody diluted at 1: 1000