

Anti-beta Actin Mouse antibody



Catalog Number: 169092

Product name

Anti-beta Actin Mouse antibody

Specificity

Human, Mouse, Monkey, Rat

Antibody description

beta Actin Mouse monoclonal antibody

Preparation

Antigen: Purified recombinant human beta actin protein fragments expressed in E.coli.

Formulation

PBS(pH 7.4) containing with 0.03% Proclin300 and 50% glycerol.

Storage

Store at 4°C short term. Store at -20°C long term. Avoid freeze / thaw cycle.

Clonality

Monoclonal

Ig Type

Mouse IgG2b


Applications

WB

Dilutions

WB: 1/10000;Other 1/5000

Validations



Western blot detection of beta actin in COS7, Jurkat, C6 and 3T3 cell lysates using beta actin mouse mAb (1:5000 diluted). Predicted band size: 45 kDa. Observed band size: 45 kDa.