

Human XEDAR/EDA2R (His Tag) recombinant protein



Catalog Number: 500492

General Information

Gene Name Synonym

X-linked ectodysplasin-A2 receptor

Protein Construction

A DNA sequence encoding the human EDA2R (NP_068555.1) extracellular domain (Met 1-Thr 138) fused with a polyhistidine tag at the C-terminus and a signal peptide at the N-terminus.

Organism

Human

Expression Host

Human Cells

QC Testing

Purity

> 95 % as determined by SDS-PAGE

Endotoxin

< 1.0 EU per μg of the protein as determined by the LAL method

Stability

Samples are stable for up to twelve months from date of receipt at -70°C

Predicted N terminal

Met

Molecular Mass

The recombinant human EDA2R consists of 149

amino acids and predicts a molecular mass of 16.7 kDa. In SDS-PAGE under reducing conditions, the apparent molecular mass of rhEDA2R is approximately 33 kDa due to glycosylation.

Formulation

Lyophilized from sterile PBS, pH 7.4

1. 5 % trehalose and mannitol are added as protectants before lyophilization.

2. Please contact us for any concerns or special requirements.

Usage Guide

Storage


Store it under sterile conditions at -20°C to -80°C .

It is recommended that the protein be aliquoted for optimal storage. Avoid repeated freeze-thaw cycles.

Reconstitution

Adding sterile water, prepare a stock solution of 0.25 mg/ml. Concentration is measured by UV-Vis.

SDS-PAGE



Human XEDAR / EDA2R Protein (His Tag) SDS-PAGE