

Anti-SRGN / serglycin antibody



Catalog Number: 102492

Product name

Anti-SRGN / serglycin antibody

Immunogen

[Human SRGN / serglycin \(His Tag\) recombinant protein](#)

Specificity

Human SRGN / serglycin

Antibody description

Rabbit polyclonal to SRGN / serglycin

Preparation

Produced in rabbits immunized with purified, recombinant Human SRGN / serglycin (rh SRGN / serglycin; CAA34900.1; Met1-Leu158). SRGN / serglycin specific IgG was purified by Human SRGN / serglycin affinity chromatography.

Formulation

0.2 µm filtered solution in PBS

Storage

This antibody can be stored at 2°C-8°C for one month without detectable loss of activity. Antibody products are stable for twelve months from date of receipt when stored at -20°C to -80°C.

Preservative-Free.

Sodium azide is recommended to avoid contamination (final concentration 0.05%-0.1%). It is toxic to cells and should be disposed of

properly. Avoid repeated freeze-thaw cycles.

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications

ELISA, WB


Dilutions

WB: 10-50 µg/ml

ELISA: 0.1-0.2 µg/mL

This antibody can be used at 0.1-0.2 µg/mL with the appropriate secondary reagents to detect Human SRGN / serglycin. The detection limit for Human SRGN / serglycin is < 0.039 ng/well.

Validations



Lanes	A
Items	A
Sample (whole cell lysate)	293T
Sample Volume (µg/lane)	30
Gel	13% SDS-PAGE reducing gel
Recommended Concentration	10-50 µg/ml
Secondary Antibody	Dylight 800 labeled Antibody to Rabbit IgG (H+L), at 1:5000 dilution.
Developed using Odyssey imaging system.	
Explanation	Predicted band size : 18 kDa Observed band size : 16 kDa

SRGN / serglycin Antibody, Rabbit PAb, Antigen Affinity Purified, Western blot