

Anti-AGBL2 antibody



Catalog Number: 175892

Product name

Polyclonal

Anti-AGBL2 antibody

Ig Type

Rabbit IgG

Specificity

Human, Mouse

Applications

WB, IHC-P

Antibody description

Rabbit polyclonal antibody to AGBL2

Dilutions

WB: 1:500-1:2000

IHC-P: 1:100-1:300

Preparation

This antigen of this antibody was the antiserum was produced against synthesized peptide derived from human cbcp2. aa range:731-780

Validations

Formulation

Liquid, Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.

Storage

Store at +4°C after thawing. Aliquot store at -20°C or -80°C. Avoid repeated freeze / thaw cycles.

Clonality




Fig1:: Western Blot analysis of various cells using CCP2 Polyclonal Antibody diluted at 1: 2000