

TBX3 antibody

Catalog Number: 115892

Product name

TBX3 antibody

Specificity

Human, Mouse, Rat; other species not tested.

Antibody description

TBX3 Rabbit Polyclonal antibody. Positive IP detected in human placenta tissue. Positive WB detected in mouse placenta tissue, HepG2 cells. Observed molecular weight by Western-blot: 80 kDa

Preparation

This antibody was obtained by immunization of TBX3 recombinant protein (Accession Number: NM_016569). Purification method: Antigen affinity purified.

Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications

ELISA, WB, IP


Dilutions

Recommended Dilution:


WB: 1:200-1:2000

IP: 1:200-1:2000

Validations



mouse placenta tissue were subjected to SDS PAGE followed by western blot with Catalog No:115892(TBX3 antibody) at dilution of 1:500



IP Result of anti-TBX3 (IP:Catalog No:115892, 4ug; Detection:Catalog No:115892 1:500) with human placenta tissue lysate 1520ug.