

Olinvacimab (Synonyms: TANIBIRUMAB)



Catalog Number: 176892

Product name

Monoclonal

Olinvacimab (Synonyms: TANIBIRUMAB)

Ig Type

Specificity

IgG1

Human VEGFR2 / KDR / CD309

Applications

Antibody description


ELISA

Anti-VEGFR2 / KDR / CD309 Antibody
(Olinvacimab)

Validations

Preparation

Recombinant expression and purified from CHO cells.



Formulation

0.1 M Pro-Ac, 20 mM Arg, pH 5.0

Storage

-80°C for 2 years under sterile conditions; -20°C for 1 year under sterile conditions; Avoid repeated freeze-thaw cycles.

SDS-PAGE

Anti-VEGFR2 / KDR / CD309 Antibody (olinvacimab) on SDS-PAGE under reducing (R) condition. The purity of the protein is greater than 95%.

Clonality