

Human ASB13 (His tag) recombinant protein



Catalog Number: 518892

General Information

Gene Name Synonym

Ankyrin repeat and SOCS box protein 13, ASB13,

Protein Construction

A DNA sequence encoding the Human ASB13 Met1-Asn278 is expressed with a 6His tag at the N-terminus.

Organism

Human

Expression Host

E. coli

QC Testing

Purity

>95% as determined by SDS-PAGE.

Endotoxin

<1.0 EU per μg protein as determined by the LAL method.

Stability

Store for up to 6 months at -20°C

Formulation


Supplied as a $0.2\ \mu\text{m}$ filtered solution of PBS, pH 7.4, 4M Urea.

Usage Guide

Storage

Store at $<-20^{\circ}\text{C}$, minimize freeze-thaw cycles.

SDS-PAGE



Human ASB13 (His tag) recombinant protein