

DEPDC6 antibody



Catalog Number: 109892

Product name

DEPDC6 antibody

IHC: 1:20-1:200

IF: 1:10-1:100


Specificity

Human, Mouse, Rat; other species not tested.

Validations

Antibody description

DEPDC6 Rabbit Polyclonal antibody. Positive WB detected in HeLa cells, A2780 cells, COLO 320 cells, HEK-293 cells, MCF7 cells. Positive IHC detected in human liver cancer. Positive IF detected in HepG2 cells. Observed molecular weight by Western-blot: 40-45kd




Preparation

This antibody was obtained by immunization of Peptide (Accession Number: NM_022783). Purification method: Antigen affinity purified.

HeLa cells were subjected to SDS PAGE followed by western blot with Catalog No:109892(DEPDC6 antibody) at dilution of 1:1000

Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.



Storage

Store at -20°C. DO NOT ALIQUOT


Clonality

Polyclonal

Immunohistochemical of paraffin-embedded human liver cancer using Catalog No:109892(DEPDC6 antibody) at dilution of 1:50 (under 10x lens)

Ig Type

Rabbit IgG



Applications

ELISA, WB, IHC, IF

Dilutions

Recommended Dilution:

WB: 1:500-1:5000


Immunohistochemical of paraffin-embedded

DEPDC6 antibody



Catalog Number: 109892

human liver cancer using Catalog
No:109892(DEPDC6 antibody) at dilution of 1:50
(under 40x lens)



Immunofluorescent analysis of HepG2 cells using
Catalog No:109892(DEPDC6 Antibody) at dilution
of 1:25 and Rhodamine-Goat anti-Rabbit IgG