

# ERAB antibody



Catalog Number: 107192

---

## Product name

ERAB antibody

## Specificity

Human, Mouse; other species not tested.

## Antibody description

ERAB Mouse Monoclonal antibody. Positive WB detected in 293 cells, HEK-293 cells, human brain tissue, human liver tissue. Observed molecular weight by Western-blot: 27kd

## Preparation

This antibody was obtained by immunization of ERAB recombinant protein (Accession Number: NM\_001037811). Purification method: Caprylic acid/ammonium sulfate precipitation.

## Formulation

PBS with 0.1% sodium azide and 50% glycerol pH 7.3.

## Storage

Store at -20°C. DO NOT ALIQUOT

## Clonality

Monoclonal

## Ig Type

Mouse IgG1

## Applications


ELISA, WB

## Dilutions

Recommended Dilution:

WB: 1:500-1:5000

## Validations



HEK-293 cells were subjected to SDS PAGE followed by western blot with Catalog No:107192(HSD17B10 antibody) at dilution of 1:1000