

# Human LILRB2/ILT4/LIR2/MIR10 (His tag) recombinant protein



Catalog Number: 519592

## General Information

### Gene Name Synonym

Leukocyte Immunoglobulin-Like Receptor Subfamily B Member 2, LIR-2, Leukocyte Immunoglobulin-Like Receptor 2, CD85 Antigen-Like Family Member D, Immunoglobulin-Like Transcript 4, ILT-4, Monocyte/Macrophage Immunoglobulin-Like Receptor 10, MIR-10, CD85d, LILRB2, ILT4, LIR2, MIR10

### Protein Construction

A DNA sequence encoding the Human LILRB2/ILT4/LIR2/MIR10 Gln22-His458 is expressed with a 6His tag at the C-terminus.

### Organism

Human

### Expression Host

Human Cells

## QC Testing

### Purity

>95% as determined by SDS-PAGE.

### Endotoxin

<1.0 EU per  $\mu\text{g}$  protein as determined by the LAL method.

### Stability

Store for up to 12 months at  $-20^{\circ}\text{C}$  to  $-80^{\circ}\text{C}$  as lyophilized powder. Aliquots of reconstituted

samples are stable at  $< -20^{\circ}\text{C}$  for 3 months.

### Formulation

Lyophilized from a  $0.2\ \mu\text{m}$  filtered solution of 20mM PB, 150mM NaCl, pH 7.2.

## Usage Guide


### Storage

Lyophilized protein should be stored at  $< -20^{\circ}\text{C}$ , though stable at room temperature for 3 weeks. Reconstituted protein solution can be stored at  $4-7^{\circ}\text{C}$  for 2-7 days.

### Reconstitution

Always centrifuge tubes before opening. Do not mix by vortex or pipetting. It is not recommended to reconstitute to a concentration less than  $100\ \mu\text{g}/\text{ml}$ . Dissolve the lyophilized protein in ddH<sub>2</sub>O. Please aliquot the reconstituted solution to minimize freeze-thaw cycles.

## SDS-PAGE



Human LILRB2/ILT4/LIR2/MIR10 (His tag)  
recombinant protein