

# MEF2C-Specific antibody



Catalog Number: 106992

## Product name

MEF2C-Specific antibody

## Specificity

Human, Mouse, Rat; other species not tested.

## Antibody description

MEF2C-Specific Rabbit Polyclonal antibody.  
Positive IHC detected in human testis tissue, human brain tissue, human kidney tissue, human lung tissue, human spleen tissue.

## Preparation

This antibody was obtained by immunization of Peptide (Accession Number: NM\_002397).  
Purification method: Antigen affinity purified.

## Formulation

PBS with 0.02% sodium azide and 50% glycerol  
pH 7.3.

## Storage

Store at -20°C. DO NOT ALIQUOT

## Clonality

Polyclonal

## Ig Type

Rabbit IgG

## Applications


ELISA, IHC

## Dilutions

Recommended Dilution:

IHC: 1:20-1:200

## Validations



Immunohistochemical of paraffin-embedded human testis using Catalog No:112637(MEF2C-Specific antibody) at dilution of 1:50 (under 10x lens)