

# Vulinacimab (Synonyms: MAB HUMAN (IGG1) ANTI P35968 (VGFR2\_HUMAN) (HLX06))



Catalog Number: 176992

## Product name

## Applications

Vulinacimab (Synonyms: MAB HUMAN (IGG1) ANTI P35968 (VGFR2\_HUMAN) (HLX06))

ELISA

## Specificity

Human VEGFR2 / KDR / CD309

## Antibody description

Anti-VEGFR2 / KDR / CD309 Antibody (Vulinacimab)

## Preparation

Recombinant expression and purified from CHO cells.

## Formulation

0.1 M Pro-Ac, 20 mM Arg, pH 5.0

## Storage

-80°C for 2 years under sterile conditions; -20°C for 1 year under sterile conditions; Avoid repeated freeze-thaw cycles.


## Clonality

Monoclonal

## Ig Type


IgG1

## Validations



## SDS-PAGE

Anti-VEGFR2 / KDR / CD309 Antibody (vulinacimab) on SDS-PAGE under reducing (R) condition. The purity of the protein is greater than 90%.



## SEC-HPLC

The purity of Anti-VEGFR2 / KDR / CD309 Antibody (vulinacimab) is more than 95%, determined by SEC-HPLC.