

DLG5 antibody

Catalog Number: 109992

Product name

DLG5 antibody

Specificity

human, mouse; other species not tested.

Antibody description

DLG5 Rabbit Polyclonal antibody. Positive WB detected in PC-3 cells, mouse kidney tissue. Observed molecular weight by Western-blot: 75kd,105kd

Preparation

This antibody was obtained by immunization of DLG5 recombinant protein (Accession Number: BC002915). Purification method: Antigen affinity purified.

Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications


ELISA, WB

Dilutions

Recommended Dilution:

WB: 1:200-1:2000

Validations



PC-3 cells were subjected to SDS PAGE followed by western blot with Catalog No:109992(DLG5 antibody) at dilution of 1:500