

Hephaestin antibody



Catalog Number: 111292

Product name

Hephaestin antibody

Specificity

Human, Mouse, Rat; other species not tested.

Antibody description

Hephaestin Rabbit Polyclonal antibody. Positive IHC detected in human pancreas cancer tissue.

Preparation

This antibody was obtained by immunization of Hephaestin recombinant protein (Accession Number: NM_001367233). Purification method: Antigen affinity purified.

Formulation

PBS with 0.1% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications


ELISA, IHC

Dilutions

Recommended Dilution:

IHC: 1:20-1:200

Validations



Immunohistochemical of paraffin-embedded human pancreas cancer using Catalog No:111292(HEPH antibody) at dilution of 1:50 (under 10x lens)