

METTL2A Polyclonal Antibody



Catalog Number: 165292

Product name

METTL2A Polyclonal Antibody

Specificity

Human

Antibody description

Polyclonal antibody to METTL2A

Preparation

Antigen: Recombinant fusion protein containing a sequence corresponding to amino acids 139-378 of human METTL2A (NP_859076.3).

Formulation

PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Storage

Store at -20°C. Avoid freeze / thaw cycles.

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications


WB, IF

Dilutions

WB 1:500 - 1:2000

IF 1:50 - 1:100

Validations



Immunofluorescence - METTL2A Polyclonal Antibody

Immunofluorescence analysis of MCF-7 cells using METTL2A antibody . Blue: DAPI for nuclear staining.