

SLAIN1 antibody



Catalog Number: 115292

Product name

SLAIN1 antibody

Specificity

Human, Mouse; other species not tested.

Antibody description

SLAIN1 Rabbit Polyclonal antibody. Positive WB detected in HeLa cells, HEK-293 cells, mouse brain tissue. Positive IF detected in HeLa cells, HEK-293 cells. Observed molecular weight by Western-blot: 65-70, 34-40 kDa

Preparation

This antibody was obtained by immunization of SLAIN1 recombinant protein (Accession Number: BC045177). Purification method: Antigen affinity purified.

Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications

ELISA, WB, IF


Dilutions

Recommended Dilution:


WB: 1:200-1:2000

IF: 1:20-1:200

Validations



HeLa cells were subjected to SDS PAGE followed by western blot with Catalog No:115292(SLAIN1 Antibody) at dilution of 1:600



Immunofluorescent analysis of HeLa cells using Catalog No:115292(SLAIN1 Antibody) at dilution of 1:50 and Alexa Fluor 488-conjugated AffiniPure Goat Anti-Rabbit IgG(H+L)