

# TNFRSF6B antibody

Catalog Number: 116292

## Product name

TNFRSF6B antibody

## Specificity

Human; other species not tested.

## Antibody description

TNFRSF6B Rabbit Polyclonal antibody. Positive IP detected in HeLa cells. Observed molecular weight by Western-blot: 32 kDa

## Preparation

This antibody was obtained by immunization of TNFRSF6B recombinant protein (Accession Number: NM\_003823). Purification method: Antigen affinity purified.

## Formulation

PBS with 0.1% sodium azide and 50% glycerol pH 7.3.

## Storage

Store at -20°C. DO NOT ALIQUOT

## Clonality

Polyclonal

## Ig Type

Rabbit IgG

## Applications


ELISA, IP

## Dilutions

Recommended Dilution:

IP: 1:200-1:1000

## Validations



IP Result of anti-TNFRSF6B (IP:Catalog No:116292, 4ug; Detection:Catalog No:116292 1:300) with HeLa cells lysate 2480ug.