

# FOXC1 antibody



Catalog Number: 110692

## Product name

FOXC1 antibody

## Specificity

Human, Mouse; other species not tested.

## Antibody description

FOXC1 Rabbit Polyclonal antibody. Positive WB detected in MCF-7 cells, HEK-293 cells, HepG2 cells. Observed molecular weight by Western-blot: 70 kDa

## Preparation

This antibody was obtained by immunization of Peptide (Accession Number: NM\_001453). Purification method: Antigen affinity purified.

## Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

## Storage

Store at -20°C. DO NOT ALIQUOT

## Clonality

Polyclonal

## Ig Type

Rabbit IgG

## Applications


ELISA, WB

## Dilutions

Recommended Dilution:

WB: 1:200-1:1000

## Validations



MCF-7 cells were subjected to SDS PAGE followed by western blot with Catalog No:110692(FOXC1 Antibody) at dilution of 1:300