

Anti-IL-4 Rabbit antibody



Catalog Number: 168692

Product name

Anti-IL-4 Rabbit antibody

Specificity

Human

Antibody description

IL-4 Rabbit monoclonal antibody

Preparation

Antigen: A synthetic peptide of human IL4

Formulation

50nM Tris-Glycine(pH 7.4), 0.15M NaCl, 40%Glycerol, 0.01% sodium azide and 0.05% BSA.

Storage

Store at 4°C short term. Store at -20°C long term. Avoid freeze / thaw cycle.

Clonality

Monoclonal

Ig Type

Rabbit IgG

Applications

WB, ICC/IF, FC


Dilutions

WB: 1/1000

ICC/IF: 1/50

FC: 1/20

Validations



Western blot detection of IL-4 in GST-IL-4 (aa25-153) cell lysates using IL-4 Rabbit mAb(1:1000 diluted). Predicted band size:17kDa. Observed band size:17kDa.