

CYP2C8 antibody



Catalog Number: 109692

Product name

CYP2C8 antibody

Specificity

Human; other species not tested.

Antibody description

CYP2C8 Rabbit Polyclonal antibody. Positive IHC detected in human liver tissue.

Preparation

This antibody was obtained by immunization of CYP2C8 recombinant protein (Accession Number: NM_000770). Purification method: Antigen affinity purified.

Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications


ELISA, IHC

Dilutions


Recommended Dilution:

IHC: 1:20-1:200

Validations



Immunohistochemical of paraffin-embedded human liver using Catalog No:109692(CYP2C8 antibody) at dilution of 1:50 (under 10x lens)



Immunohistochemical of paraffin-embedded human liver using Catalog No:109692(CYP2C8 antibody) at dilution of 1:50 (under 40x lens)