

Human TINAGL1/GIS5/LCN7/OLRG2 (His tag) recombinant protein



Catalog Number: 519692

General Information

Gene Name Synonym

Tubulointerstitial nephritis antigen-like (TINAGL1), also known as Glucocorticoid-inducible protein 5, Oxidized LDL-responsive gene 2 protein, TIN Ag-related protein, Tubulointerstitial nephritis antigen-related protein, GIS5, LCN7, OLRG2 and TINAGL. belongs to the peptidase C1 family.

Protein Construction

A DNA sequence encoding the Human TINAGL1/GIS5/LCN7/OLRG2 Ala22-His467 is expressed with a 6His tag at the C-terminus.

Organism

Human

Expression Host

Human Cells

QC Testing

Purity

>95% as determined by SDS-PAGE.

Endotoxin

<1.0 EU per μg protein as determined by the LAL method.

Stability

Store for up to 6 months at -20°C

Formulation


Supplied as a $0.2\ \mu\text{m}$ filtered solution of 20mM PB, 150mM NaCl, 1mM EDTA, 2m MDTT, pH7.4.

Usage Guide

Storage

Store at $<-20^{\circ}\text{C}$, minimize freeze-thaw cycles.

SDS-PAGE



Human TINAGL1/GIS5/LCN7/OLRG2 (His tag)
recombinant protein