

PSMC6 antibody

Catalog Number: 114392

Product name

PSMC6 antibody

Specificity

Human, Mouse, Rat; other species not tested.

Antibody description

PSMC6 Rabbit Polyclonal antibody. Positive WB detected in human brain tissue, human testis tissue, mouse brain tissue, rat liver tissue. Positive IF detected in SH-SY5Y cells. Observed molecular weight by Western-blot: 42 kDa

Preparation

This antibody was obtained by immunization of PSMC6 recombinant protein (Accession Number: NM_002806). Purification method: Antigen affinity purified.

Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications

ELISA, WB, IF


Dilutions

Recommended Dilution:


WB: 1:1000-1:10000

IF: 1:20-1:200

Validations



human brain tissue were subjected to SDS PAGE followed by western blot with Catalog No:114392(PSMC6 antibody) at dilution of 1:500



Immunofluorescent analysis of (-20oc Ethanol) fixed SH-SY5Y cells using Catalog No:114392(PSMC6 Antibody) at dilution of 1:50 and Alexa Fluor 488-conjugated AffiniPure Goat Anti-Rabbit IgG(H+L)