

IL21R antibody

Catalog Number: 107392

Product name

IL21R antibody

Specificity

Human; other species not tested.

Antibody description

IL21R Mouse Monoclonal antibody. Positive IHC detected in human tonsillitis tissue. Positive WB detected in K-562 cells, Jurkat cells. Observed molecular weight by Western-blot: 65 kDa

Preparation

This antibody was obtained by immunization of IL21R recombinant protein (Accession Number: NM_181078). Purification method: Protein A purified.

Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Monoclonal

Ig Type

Mouse IgG2b

Applications

ELISA, IHC, WB


Dilutions

Recommended Dilution:


WB: 1:500-1:5000

IHC: 1:50-1:200


Validations



K-562 cells were subjected to SDS PAGE followed by western blot with Catalog No:107392(IL21R Antibody) at dilution of 1:1500



Immunohistochemistry of paraffin-embedded human tonsillitis tissue slide using (IL21R Antibody) at dilution of 1:100 (under 10x lens)



Immunohistochemistry of paraffin-embedded human tonsillitis tissue slide using (IL21R Antibody) at dilution of 1:100 (under 40x lens)