

Anti-SKP2 (4D8) Mouse antibody



Catalog Number: 168392

Product name

Anti-SKP2 (4D8) Mouse antibody

Specificity

Human

Antibody description

SKP2 (4D8) Mouse monoclonal antibody

Preparation

Antigen: Purified recombinant fragment of SKP2 (aa1-130) expressed in E. Coli.

Formulation

Ascitic fluid containing 0.03% sodium azide.

Storage

Store at 4°C short term. Store at -20°C long term. Avoid freeze / thaw cycle.

Clonality

Monoclonal

Ig Type

Mouse IgG1

Applications


WB

Dilutions

WB: 1/500 - 1/2000

ELISA: 1/10000

Validations



Western blot detection of Skp2 in HeLa cell lysates using Skp2 mouse mAb (1:2000 diluted). Predicted band size: 48KDa. Observed band size: 48KDa.