

B4GALNT2 antibody



Catalog Number: 108392

Product name

B4GALNT2 antibody

Specificity

Human; other species not tested.

Antibody description

B4GALNT2 Rabbit Polyclonal antibody. Positive WB detected in Jurkat cells, COLO 320 cells. Positive IHC detected in human colon cancer tissue, human colon tissue. Observed molecular weight by Western-blot: 65 kDa

Preparation

This antibody was obtained by immunization of B4GALNT2 recombinant protein (Accession Number: NM_153446). Purification method: Antigen Affinity purified.

Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications

ELISA, WB, IHC


Dilutions

Recommended Dilution:


WB: 1:500-1:5000

IHC: 1:20-1:200


Validations



Jurkat cells were subjected to SDS PAGE followed by western blot with Catalog No:108392(B4GALNT2 Antibody) at dilution of



Immunohistochemistry of paraffin-embedded human colon cancer tissue slide using Catalog No:108392(B4GALNT2 Antibody) at dilution of 1:50 (under 10x lens)



Immunohistochemistry of paraffin-embedded human colon cancer tissue slide using Catalog

B4GALNT2 antibody



Catalog Number: 108392

No:108392(B4GALNT2 Antibody) at dilution of 1:50 (under 40x lens)