

# FUK antibody



Catalog Number: 110792

## Product name

FUK antibody

## Specificity

Human, Mouse, Rat; other species not tested.

## Antibody description

FUK Rabbit Polyclonal antibody. Positive IHC detected in mouse kidney tissue. Positive WB detected in mouse kidney tissue, mouse thymus tissue. Observed molecular weight by Western-blot: 107kd

## Preparation

This antibody was obtained by immunization of FUK recombinant protein (Accession Number: XM\_011522929). Purification method: Antigen affinity purified.

## Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

## Storage

Store at -20°C. DO NOT ALIQUOT

## Clonality

Polyclonal

## Ig Type

Rabbit IgG

## Applications

ELISA, WB, IHC


## Dilutions

Recommended Dilution:


WB: 1:500-1:5000

IHC: 1:20-1:200

## Validations



mouse kidney tissue were subjected to SDS PAGE followed by western blot with Catalog No:110792(FUK antibody) at dilution of 1:500



Immunohistochemistry of paraffin-embedded mouse kidney tissue slide using Catalog No:110792(FUK Antibody) at dilution of 1:50 (under 40x lens)