

Mouse CD33 / Siglec-3 (Fc Tag) recombinant protein



Catalog Number: 502792

General Information

Gene Name Synonym

Sialic acid-binding Ig-like lectin 3

Protein Construction

A DNA sequence encoding the extracellular domain of mouse CD33 (Q63994-1) (Met 1-Glu 240) was fused with the Fc region of human IgG1 at the C-terminus.

Organism

Mouse

Expression Host

Human Cells

QC Testing

Purity

> 90 % as determined by SDS-PAGE

Endotoxin

< 1.0 EU per μg of the protein as determined by the LAL method

Stability

Samples are stable for up to twelve months from date of receipt at -70°C

Predicted N terminal

Gln 17

Molecular Mass

The secreted recombinant mouse CD33/Fc is a disulfide-linked homodimer. The reduced

monomer comprises 465 amino acids and has a calculated molecular mass of 51.6 kDa. As a result of glycosylation, the apparent molecular mass of rm CD33/Fc monomer is approximately 55-60 kDa in SDS-PAGE under reducing conditions.

Formulation

Lyophilized from sterile PBS, pH 7.4

1. 5 % trehalose and mannitol are added as protectants before lyophilization.
2. Please contact us for any concerns or special requirements.

Usage Guide


Storage

Store it under sterile conditions at -20°C to -80°C . It is recommended that the protein be aliquoted for optimal storage. Avoid repeated freeze-thaw cycles.

Reconstitution

Adding sterile water, prepare a stock solution of 0.25 mg/ml. Concentration is measured by UV-Vis.

SDS-PAGE



Mouse CD33 / Siglec-3 Protein (Fc Tag) SDS-PAGE