

Human ACE2 Protein (Fc Tag)



Catalog Number: 521003

General Information

Gene Name Synonym

Angiotensin-Converting Enzyme 2

Protein Construction

Fc tag was fused to the C-terminus of Human ACE2 Protein

Organism

Human

Expression Host

Human Cells

QC Testing

Purity

>90%

Endotoxin

<5 EU/mg

Stability

Samples are stable for up to twelve months from date of receipt at -70°C

Molecular Mass

110 kDa (The protein migrates as a higher molecular weight under reducing (R) condition (SDS-PAGE) due to glycosylation)

Formulation


Liquid, buffer: PBS(pH7.2)

Usage Guide

Storage

Store it under sterile conditions at -20°C to -80°C. It is recommended that the protein be aliquoted for optimal storage. Avoid repeated freeze-thaw cycles.

SDS-PAGE



Human ACE2 Protein (Fc Tag)