

GNGT1 antibody



Catalog Number: 111003

Product name

GNGT1 antibody

Specificity

Human, Mouse, Rat; other species not tested.

Antibody description

GNGT1 Rabbit Polyclonal antibody. Positive WB detected in mouse eye tissue. Positive IF detected in HeLa cells. Observed molecular weight by Western-blot: 8 kDa

Preparation

This antibody was obtained by immunization of GNGT1 recombinant protein (Accession Number: NM_001329426). Purification method: Antigen affinity purified.

Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications

ELISA, WB, IF


Dilutions

Recommended Dilution:


WB: 1:200-1:1000

IF: 1:10-1:100

Validations



mouse eye tissue were subjected to SDS PAGE followed by western blot with Catalog No:111003(GNGT1 antibody) at dilution of 1:100



Immunofluorescent analysis of HeLa cells using Catalog No:111003(GNGT1 Antibody) at dilution of 1:25 and Alexa Fluor 488-conjugated AffiniPure Goat Anti-Rabbit IgG(H+L)