

# Anti-Dppa2 antibody



Catalog Number: 175003

## Product name

Anti-Dppa2 antibody

## Ig Type

Rabbit IgG

## Specificity

Mouse

## Applications

WB

## Antibody description

Rabbit polyclonal antibody to Dppa2

## Dilutions

WB: 1:500

## Preparation

This antigen of this antibody was peptide

## Validations

## Formulation

Liquid, 1\*PBS (pH7.4), 0.2% BSA, 40% Glycerol.  
Preservative: 0.05% Sodium Azide.

## Storage

Store at +4°C after thawing. Aliquot store at -20°C  
or -80°C. Avoid repeated freeze / thaw cycles.

## Clonality

Polyclonal




Fig1.; Western blot analysis on F9 using anti-DPPA2 polyclonal antibody