

# BNIP3 Polyclonal Antibody



Catalog Number: 162403

---

## Product name

BNIP3 Polyclonal Antibody

## Specificity

Human

## Antibody description

Polyclonal antibody to BNIP3

## Preparation

Antigen: Recombinant fusion protein containing a sequence corresponding to amino acids 1-140 of human BNIP3 (NP\_004043.3).

## Formulation

PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

## Storage

Store at -20°C. Avoid freeze / thaw cycles.

## Clonality

Polyclonal

## Ig Type

Rabbit IgG


## Applications

IF

## Dilutions

IF 1:50 - 1:200

## Validations



Immunofluorescence - BNIP3 Polyclonal Antibody

Immunofluorescence analysis of U2OS cells using BNIP3 antibody at dilution of 1:100. Blue: DAPI for nuclear staining.