

# LY6K antibody

Catalog Number: 112403

## Product name

LY6K antibody

## Specificity

Human; other species not tested.

## Antibody description

LY6K Rabbit Polyclonal antibody. Positive IP detected in PC-3 cells. Positive WB detected in human testis tissue. Observed molecular weight by Western-blot: 28 kDa

## Preparation

This antibody was obtained by immunization of LY6K recombinant protein (Accession Number: BC001291). Purification method: Antigen affinity purified.

## Formulation

PBS with 0.1% sodium azide and 50% glycerol pH 7.3.

## Storage

Store at -20°C. DO NOT ALIQUOT

## Clonality

Polyclonal

## Ig Type

Rabbit IgG

## Applications

ELISA, WB, IP


## Dilutions

Recommended Dilution:


WB: 1:200-1:2000

IP: 1:200-1:1000

## Validations



human testis tissue were subjected to SDS PAGE followed by western blot with Catalog No:112403(LY6K antibody) at dilution of 1:500



IP Result of anti-LY6K (IP:Catalog No:112403, 4ug; Detection:Catalog No:112403 1:300) with PC-3 cells lysate 3200ug.