

Anti-NUP88 antibody



Catalog Number: 176403

Product name

Anti-NUP88 antibody

Specificity

Human, Mouse

Antibody description

Rabbit polyclonal antibody to NUP88

Preparation

This antigen of this antibody was synthesized peptide derived from nup88 . at aa range: 330-410

Formulation

Liquid, Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.

Storage

Store at +4°C after thawing. Aliquot store at -20°C or -80°C. Avoid repeated freeze / thaw cycles.

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications

WB

Dilutions

WB: 1:500-1:2000

Validations




Fig1:; Western Blot analysis of extracts from K562 cells, using Nup88 Polyclonal Antibody.. Secondary antibody (catalog#: HA1001) was diluted at 1:20000