

Anti-HDAC6 (3B2) Mouse antibody



Catalog Number: 168403

Product name

Anti-HDAC6 (3B2) Mouse antibody

Ig Type

Mouse IgG2a

Specificity

Human

Applications

WB

Antibody description

HDAC6 (3B2) Mouse monoclonal antibody

Dilutions

WB: 1/1000

Preparation

Antigen: Purified recombinant human HDAC6 protein expressed in E.coli.

Formulation

Purified mouse monoclonal antibody in PBS(pH 7.4) containing with 0.03% Proclin300 and 50% glycerol.


Storage

Store at 4°C short term. Store at -20°C long term. Avoid freeze / thaw cycle.

Clonality

Monoclonal

Validations



Western blot detection of HDAC6 in T47D, Jurkat CHO-K1(B) and CHO-K1 transfected by HDAC6-fragment EGFP fusion proteinuff08Auff09cell lysates using HDAC6 mouse mAb (1:1000 diluted).