

# Anti-HDAC6 (3B2) Mouse antibody



Catalog Number: 168403

## Product name

Anti-HDAC6 (3B2) Mouse antibody

## Specificity

Human

## Antibody description

HDAC6 (3B2) Mouse monoclonal antibody

## Preparation

Antigen: Purified recombinant human HDAC6 protein expressed in E.coli.

## Formulation

Purified mouse monoclonal antibody in PBS(pH 7.4) containing with 0.03% Proclin300 and 50% glycerol.

## Storage

Store at 4°C short term. Store at -20°C long term. Avoid freeze / thaw cycle.

## Clonality

Monoclonal

## Ig Type

Mouse IgG2a


## Applications

WB

## Dilutions

WB: 1/1000

## Validations



Western blot detection of HDAC6 in T47D, Jurkat CHO-K1(B) and CHO-K1 transfected by HDAC6-fragment EGFP fusion proteinuff08Auff09cell lysates using HDAC6 mouse mAb (1:1000 diluted).