

# ABCD2 antibody

Catalog Number: 107803

## Product name

ABCD2 antibody

## Specificity

Human, Mouse, Rat; other species not tested.

## Antibody description

ABCD2 Rabbit Polyclonal antibody. Positive WB detected in mouse cerebellum tissue. Positive IP detected in A431 cells. Observed molecular weight by Western-blot: 75 kDa

## Preparation

This antibody was obtained by immunization of ABCD2 recombinant protein (Accession Number: NM\_005164). Purification method: Antigen affinity purified.

## Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

## Storage

Store at -20°C. DO NOT ALIQUOT

## Clonality

Polyclonal

## Ig Type

Rabbit IgG

## Applications

ELISA, WB, IP


## Dilutions

Recommended Dilution:


WB: 1:200-1:1000

IP: 1:200-1:1000

## Validations



mouse cerebellum tissue were subjected to SDS PAGE followed by western blot with Catalog No:107803(ABCD2 Antibody) at dilution of 1:300



IP Result of anti-ABCD2 (IP:Catalog No:107803, 4ug; Detection:Catalog No:107803 1:300) with A431 cells lysate 1200ug.