

LIMK1 Monoclonal Antibody



Catalog Number: 160103

Product name

LIMK1 Monoclonal Antibody

Specificity

Human, Mouse

Antibody description

Monoclonal antibody to LIMK1

Preparation

Antigen: A synthetic peptide of human LIMK1

Formulation

PBS with 0.02% sodium azide, pH7.3.

Storage

Store at 4°C. Avoid freeze / thaw cycles.

Clonality

Monoclonal

Ig Type

Mouse IgG


Applications

WB

Dilutions

WB 1:500 - 1:1000

Validations



Western blot - LIMK1 Monoclonal Antibody

Western blot analysis of extracts of HeLa cells, using LIMK1 antibody. Secondary antibody: HRP Goat Anti-Mouse IgG (H+L) at 1:10000 dilution. Lysates/proteins: 25ug per lane. Blocking buffer: 3% nonfat dry milk in TBST.