

Anti-CES3/Carboxylesterase 3 antibody



Catalog Number: 100103

Product name

Anti-CES3/Carboxylesterase 3 antibody

Specificity

Human CES3 / Carboxylesterase 3

Antibody description

Rabbit monoclonal to CES3/Carboxylesterase 3

Preparation

This antibody was obtained from a rabbit immunized with purified, recombinant Human CES3 / Carboxylesterase 3 (rh CES3 / Carboxylesterase 3; NP_079198.2; Met1-Glu561).

Formulation

0.2 µm filtered solution in PBS

Storage

This antibody can be stored at 2°C-8°C for one month without detectable loss of activity. Antibody products are stable for twelve months from date of receipt when stored at -20°C to -80°C.

Preservative-Free.

Sodium azide is recommended to avoid contamination (final concentration 0.05%-0.1%). It is toxic to cells and should be disposed of properly. Avoid repeated freeze-thaw cycles.

Clonality

Monoclonal

Ig Type

Rabbit IgG


Applications

IHC-P

Dilutions


IHC-P: 1-10 µg/mL

Validations



CES3 / Carboxylesterase 3 Antibody, Rabbit MAb, Immunohistochemistry

Immunochemical staining of human CES3 in human liver with rabbit monoclonal antibody (1 µg/mL, formalin-fixed paraffin embedded



sections).

CES3 / Carboxylesterase 3 Antibody, Rabbit MAb, Immunohistochemistry

Immunochemical staining of human CES3 in human colon with rabbit monoclonal antibody (1 µg/mL, formalin-fixed paraffin embedded sections).