

GMPR antibody

Catalog Number: 111103

Product name

GMPR antibody

Specificity

human, mouse, rat; other species not tested.

Antibody description

GMPR Rabbit Polyclonal antibody. Positive IHC detected in human colon cancer. Positive WB detected in mouse skin tissue, mouse skeletal muscle tissue. Observed molecular weight by Western-blot: 37kd

Preparation

This antibody was obtained by immunization of GMPR recombinant protein (Accession Number: NM_006877). Purification method: Antigen affinity purified.

Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications

ELISA, WB, IHC


Dilutions

Recommended Dilution:


WB: 1:200-1:2000

IHC: 1:20-1:200


Validations



mouse skin tissue were subjected to SDS PAGE followed by western blot with Catalog No:111103(GMPR antibody) at dilution of 1:500



Immunohistochemical of paraffin-embedded human colon cancer using Catalog No:111103(GMPR antibody) at dilution of 1:100 (under 10x lens)



Immunohistochemical of paraffin-embedded human colon cancer using Catalog No:111103(GMPR antibody) at dilution of 1:100

GMPR antibody



Catalog Number: 111103

(under 40x lens)