

NEK2 antibody



Catalog Number: 113103

Product name

NEK2 antibody

Specificity

human, mouse, rat; other species not tested.

Antibody description

NEK2 Rabbit Polyclonal antibody. Positive WB detected in K-562 cells. Positive IF detected in Hela. Observed molecular weight by Western-blot: 48kd

Preparation

This antibody was obtained by immunization of NEK2 recombinant protein (Accession Number: NM_002497). Purification method: Antigen affinity purified.

Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications

ELISA, WB, IF


Dilutions

Recommended Dilution:


WB: 1:200-1:2000

IF: 1:10-1:100

Validations



K-562 cells were subjected to SDS PAGE followed by western blot with Catalog No:113103(NEK2 antibody) at dilution of 1:500



Immunofluorescent analysis of Hela cells, using NEK2 antibody Catalog No:113103 at 1:25 dilution and Rhodamine-labeled goat anti-rabbit IgG (red).