

PPP2R3B antibody



Catalog Number: 114103

Product name

PPP2R3B antibody

Specificity

Human; other species not tested.

Antibody description

PPP2R3B Rabbit Polyclonal antibody. Positive WB detected in HL-60 cells, HeLa cells. Observed molecular weight by Western-blot: 65 kDa

Preparation

This antibody was obtained by immunization of PPP2R3B recombinant protein (Accession Number: NM_199326). Purification method: Antigen affinity purified.

Formulation

PBS with 0.1% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications


ELISA, WB

Dilutions

Recommended Dilution:

WB: 1:500-1:5000

Validations



HL-60 cells were subjected to SDS PAGE followed by western blot with Catalog No:114103(PPP2R3B Antibody) at dilution of 1:1000