

SEMA3C antibody



Catalog Number: 115103

Product name

SEMA3C antibody

Specificity

Human, Mouse; other species not tested.

Antibody description

SEMA3C Rabbit Polyclonal antibody. Positive WB detected in mouse heart tissue. Observed molecular weight by Western-blot: 70 kDa

Preparation

This antibody was obtained by immunization of SEMA3C recombinant protein (Accession Number: NM_006379). Purification method: Antigen affinity purified.

Formulation

PBS with 0.1% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications


ELISA, WB

Dilutions

Recommended Dilution:

WB: 1:500-1:5000

Validations



mouse heart tissue were subjected to SDS PAGE followed by western blot with Catalog No:115103(SEMA3C antibody) at dilution of 1:500